

High-Speed CO₂ Ramping Technique


Rapid A/C_i Curves in Minutes

CIRAS-3 is capable of rapidly controlling CO₂ gas concentration while simultaneously and continuously recording data. This capability has been available using the multiple level response scripts in every CIRAS-3. With new interest in a linear ramp of CO₂ concentration, PP Systems has released a new version of the PC-based Script Editor that makes it simpler to create the response script file needed to create linear ramp response curves, where the CO₂ reference level is changed at each recording interval.

This application note describes the straightforward steps needed to create the script file, set up the CIRAS-3, run a linear ramp experiment, and record data. Data is plotted to illustrate the linear ramp capability. Post processing of the data to generate A vs. C_i curves from the ramped gas exchange data is then described.

Script File Generation

1. Start the CIRAS-3 Utility program (**Start > All Programs > PP Systems > CIRAS-3 Utility**). Requires Version 2.00 or later to get the ramp features described in this application note.
2. Choose **File > Open Scripts Folder**. Choose a location where the response scripts are saved. The default is **User\MyDocuments\PP Systems\Ciras3_PC_Utility\Response**.
3. The initial screen that opens is the standard script editor window that allows discrete increments of CO₂, H₂O, Temperature and PAR. Click the **Ramp** checkbox to open the Ramp Editor. One ramping script for C₃ leaves, and one for C₄ leaves are available from PP Systems. The two scripts differ primarily in the lower maximum CO₂ level reached for C₄ photosynthesis. Both may be edited by the user.
4. Edit **CO₂ Start, CO₂ Stop, Ramp Time** and other experimental parameters in any of the white text boxes.



5. The **Update Levels** button will highlight in blue on any change to indicate that the grid needs updating. Click the blue **Update Levels** button to recompute the grid from current parameters.
6. The column labeled **Record Interval** shows as 1, indicating that 1 data point is recorded for each recording interval, or one data point every 1.6 seconds.
7. Select **File > Save As** to save the new response script. The name of all ramp scripts should include the word "Ramp" somewhere in the filename to allow future editing in the Ramp Editor (although any name can be saved and successfully used in the CIRAS-3).

A screenshot of the 'Script file: Default_RAMP' window. The 'Ramp' checkbox is checked. Parameters are set as follows: Number of Levels: 226; Start Acclimation: 120 (s); CO₂ Start: 200; CO₂ Stop: 1000 (umol mol⁻¹); Ramp Time: 6 (minutes); H₂O: 100 (%); Leaf Temperature: 25 (C); Light Intensity: 1000 (umol m⁻² s⁻¹); RGBW: Red 38, Green 37, Blue 25, White 0 (%). The 'Update Levels' button is highlighted in blue. A table at the bottom shows the script levels.

Level	Acclimation	Records	Record Interval	CO ₂	H ₂ O	Leaf Temp	PAR	RGBW
1	120	1	1	200	100	25	1000	38-37-25-0
2	1	1	1	204	100	25	1000	38-37-25-0
3	1	1	1	207	100	25	1000	38-37-25-0
4	1	1	1	211	100	25	1000	38-37-25-0
5	1	1	1	214	100	25	1000	38-37-25-0
6	1	1	1	218	100	25	1000	38-37-25-0
7	1	1	1	224	100	25	1000	38-37-25-0

Set Up CIRAS-3 and Record Data

Once the script file is created in the PC Utility package, follow these steps to prepare the CIRAS-3 for a linear ramp process.

1. Transfer the script file from the PC to the CIRAS-3 using a USB flash drive: Navigate to **Operations > Rec Options > Edit Rsp Crv > Transfer**, and follow prompts to Import the script file into **Internal Memory Files**.
2. Perform a **Stored Diff Bal Calibration** to allow the CIRAS-3 to have accurate offset information applied continuously throughout the linear ramp: Navigate to **Operations > Calibration > Stored Diff Bal** and follow the prompts. Set the CO₂ min and CO₂ max to match the ramp min and max concentrations. The H₂O range can be left at the default of 0 to 100%, or reduced to a smaller range closer to the H₂O operating point of the response script. Choose **Start** and the CIRAS-3 steps through 6 levels between the min and max settings and performs a diff-balance at each level, then computes a regression analysis of the resulting offsets that will then be applied for every CO₂ concentration in real time. The Stored Diff Balance Calibration takes 20-30 minutes and must be allowed to continue to completion. Under **Settings** make sure **Zero, Diff Bal Mode** is set to **Auto Zero, stored Diff Bal**.
3. Start the experiment: The first ramp recording creates the baseline trace to characterize the time response of the system and stores the data for subsequent post processing. The PLC3 Universal Leaf Cuvette should be set up with appropriate operational settings similar to the start of the ramp and the cuvette should be closed with no leaf. Navigate to **Operations > Rec Options**. Choose **Response Curves** in top line to enable the Response Curve Scripts dropdown box. Choose the correct response script for the linear ramp. Choose either internal storage or external storage (USB flash drive) for the results data files. Choose **Start** to begin the experiment. Choose **Back** to return to the data or graph screen and watch the ramp script occur. Red triangles are shown on the graph at each recorded data point (and at this recording interval, the red triangles overlap).
4. Start the experiment with a leaf: Follow the same steps as step 3, except with a leaf in the PLC3 chamber. Note, if the current environment of the leaf to be tested differs greatly from the first step in the ramping script, it is recommended to allow the leaf to equilibrate to those cuvette conditions prior to starting the ramping script with the leaf. If one watches the plot of A vs. time for the ramp with the leaf, it becomes clear when CO₂ saturation occurs, and the ramp can be terminated at that time, by hitting the **End Recording** button.
5. Transfer data files to PC and begin post processing. If data files were stored to the internal memory, use the **Operations > Rec Options > Transfer Data** screen to move the files to a USB flash drive.
6. Post Processing: Because of the Stored Diff Balance capability of the CIRAS-3, post processing can be very simple – all done in Excel by adding just 3 new columns to the standard CIRAS-3 output file. To begin, open the file for the empty cuvette ramp in Excel. Copy the column of assimilation rates (A), which is column Y. Open the file with the ramp with the leaf in the cuvette, and paste the empty cuvette A rates into an empty column. Subtract the empty cuvette A rates from the rates obtained with the leaf in the cuvette. The differences are the actual assimilation rates, starting after a brief lag period (about line 20, usually).
7. Compute C_i: The sub-stomatal CO₂ concentration, C_i, is now recomputed using the new actual Assimilation and the other values that have not changed (CO₂ analysis, and g_s and E.):


$$C_i (\mu\text{mol mol}^{-1}) = \frac{\left[\left(g_c - \frac{E}{2} \right) \times C_{\text{out}} \right] - A}{\left(g_c + \frac{E}{2} \right)}$$

Where g_c is the total conductance to CO₂ transfer:

$$g_c (\text{mmol m}^{-2} \text{ s}^{-1}) = \left[\frac{1}{(1.585 \times r_a) + (1.37 \times r_b)} \right] \times 10^3$$

[1.585 is the diffusion ratio of CO₂ and water in *air*, and 1.37 is the diffusion ratio of CO₂ and water in the *boundary layer*.]

An Excel spreadsheet is available from PP Systems to use as a template for this calculation.
8. Plot A vs. C_i, starting at about line 20, after the linear ramp stabilizes.




Notes:

1. The Stored Diff Balance Calibration should be performed at the beginning of the each day after a 30-minute initial warmup. The Stored Diff Balance is then stable throughout an 8-12 hour day of testing and does not need to be redone until the next day.
2. The No Leaf (empty chamber) ramp is stable for 4 hours or more and does not need to be redone unless experimental parameters are changed.
3. We have not encountered a situation where humidity values nor g_s values change rapidly enough during the ramp to cause substantial errors in g_s and calculated C_i. The scripts provided allow curves up to saturating A to be completed in about 5 minutes, and g_s changes with C_i are normally fairly small in that time.

Sample Results

A linear ramp response curve from 200 ppm to 600 ppm over 6 minutes was created with the Script Editor with a 300 ml/min cuvette flow, a 100 ml/min analysis flow and PAR set to 1000 $\mu\text{mol mol}^{-2} \text{s}^{-1}$. The CIRAS-3 was warmed up and a Stored Diff Balance Calibration was performed. The response script was run once with an empty chamber, and again with an *ocimum basilicum* leaf in the PLC3 chamber.




In the no-leaf case, the CO_2 differential is a relatively constant -14 ppm for most of the ramp, after starting out at 0 ppm during the 120 second acclimation time when the $\text{CO}_{2,r}$ is kept constant at 200 ppm. The -14 ppm $\text{CO}_{2,d}$ represents the response time of the system including cuvette mixing and gas transport back to the CIRAS-3 console, equivalent to 12.6 sec with these particular settings. Faster response time could be obtained with higher cuvette flow rate, but with a corresponding lower CO_2 differential. Since the CIRAS-3 can perform a stored diff-balance over the full range of the ramp prior to running the response script, the reference and analysis do not have to be corrected for accumulated channel difference.

With an active leaf in the PLC3, the $\text{CO}_{2,d}$ begins at -17 ppm during the 120 sec acclimation (instead of 0 as in the no-leaf case) because the leaf is actively assimilating. Then as the $\text{CO}_{2,r}$ increases from 200 to 600, the $\text{CO}_{2,d}$ increases from -16 ppm to -65 ppm at the end of the ramp.

A vs C_i Comparisons


The Rapid A-C_i curve technique and traditional point-by-point steady state A-C_i technique were compared on identical leaves a few minutes apart. Data on both C₃ (soybean) and a C₄ (giant foxtail) were made and show very good agreement between the two methods.

Giant Foxtail at 25°C



Comparison of High-Speed A/C_i Ramping (black points) to traditional point-by-point Steady State (red points) for a typical C₄ Giant Foxtail leaf with PAR of 1500 μmol m⁻² s⁻¹ and Cuvette Flow of 300 ml/min. Reference CO₂ was ramped from 50 to 500 in 5 minutes (with one initial 2 minute acclimation). Each Steady State point had a 2 minute acclimation time for total data recording time of 18 minutes.

Soybean at 25°C



Comparison of High-Speed A/C_i Ramping (black points) to traditional point-by-point Steady State (red points) for a typical C₃ Soybean leaf with PAR of 1500 μmol m⁻² s⁻¹ and Cuvette Flow of 300 ml/min. Reference CO₂ ramped from 100 to 1000 in 8 minutes (with one initial 2 minute acclimation). Each Steady State point had a 2 minute acclimation time for total data recording time of 22 minutes.



110 Haverhill Road, Suite 301
Amesbury, MA 01913 U.S.A.

Tel: +1 978-834-0505
Fax: +1 978-834-0545

support@ppsystems.com
ppsystems.com